

VK Solutions and Tablets®

Powerful cleaning with minimal effort.

- Powerful ultrasonic energy quickly removes and breaks down hard to clean deposits like blood, bone, and food particles that chemical detergents can leave behind.
- Powerful, yet gentle on your instruments.
- Biodegradable and pH neutral*.
- Non-toxic, non-corrosive formula contains no chlorine or harsh alkaline to damage your instruments*.

KULZER
MITSUI CHEMICALS GROUP

Giving a hand to oral health.

^{*}Refers to VibraKleen E2 Ultrasonic Cleaning Tablets.

Product	Features	Benefits	Uses	Packaging	Product Numbers
VIBRAKLEE	EN E ^{2®} ULTRASONIC CLEAI	NING SOLUTIONS			
VibraKleen E ² Ultrasonic Cleaning Solutions	 Effervescent, enzymatic cleaning tablets Non-corrosive and low-foaming Biodegradable, non-toxic & pH neutral No chlorine or harsh alkaline Drop directly into ultrasonic cleaner Compact, convenient package of 20 tablets dilutes to 10 gallons Clean mint scent Effective as evacuation cleaner Power of protease enzymes 	Breaks down hard-to-clean deposits in ultrasonic and evacuation systems Easy on instruments Environmentally friendly and safe to work with Won't damage stainless steel equipment No messy mixing of liquid concentrates Saves on counter and storage space Minimizes any odor from tank, pleasant to work with One product for multiple uses Effective, fast removal of blood, bone and food particles that chemical detergents leave behind	• For general purpose ultrasonic cleaning, evacuation system cleaning and pre-soak	Box of 20 Box of 80 2 tablets = 1 gallon	50037100— Box of 20 (Case Lot 24) 50037101— Clinic Pack— Box of 80 (Case Lot 12)
DRI-CLAVE® VK ULTRASONIC CLEANING SOLUTIONS					
VK-1 General Purpose Cleaner	 Efficient cleansing time (2–10) minutes 10:1 concentrated formulation Color coded solution Anticorrosive 	 Protects your equipment Avoid messy cleaning Economical gallon container Mistake-free identification 	Removes debris from instruments, burs, general equipment, ortho appliances, gold and porcelain Use directly in ultrasonic tank	Gallon— blue Liquid	50036202* 40000239* (Canadian customers only)
VK-3 Plaster and Stone Remover	 Efficient cleansing time (5–10) minutes No diluting necessary—use at full strength Color coded solution 	 Protects your instruments Avoid messy cleaning Saves time Economical gallon container Mistake-free identification 	 For removal of plaster and stone Do not use with aluminum Place in glass breaker inside ultra- sonic cleaner 	Gallon— pink Liquid	50036206* 40000240* (Canadian customers only)
VK-4 Tartar and Stain Remover	 Efficient cleansing time (5–10) minutes No diluting necessary— use at full strength 	 Avoid messy cleaning Ready-to-use formula Saves time Economical gallon container 	Removal of tartar and stains 5-10 minutes only Place in glass beaker inside ultrasonic cleaner	Gallon— clear Liquid	50036208* 40000241* (Canadian customers only)
VK-6 Temporary Cement Remover	 Efficient cleansing time (5-10) minutes No diluting necessary— use at full strength Color coded solution 	 Protects your instruments Saves time Economical gallon container 	Removal of temporary cement including cavit and eugenol-based temporary cement 1 – 10 minutes only Place in glass beaker inside ultrasonic cleaner	Gallon— yellow Liquid	50036212* 40000242* (Canadian customers only)
VACUKLEE	N E ^{2®} EVACUATION SYSTE	M CLEANER			
VacuKleen E ² Evacuation System Cleaner	Highly effective powdered formula Non-corrosive and low-foaming Biodegradable & non-toxic Convenient, easy-mix packets Power of protease enzymes Fresh, clean, scent of lemon Simple to use	Removes hard-to-clean protein-based deposits that soapy solutions can leave behind Environmentally friendly and safe to work with Pre-measured packets—saves time and minimizes waste Effective and quick dissolving of blood, tissue, saliva and encrusted residue in evacuation lines Real lemon extracts gives a clean smell that eliminates evacuation system odor and makes the operatories smell fresh Add powder packet to warm water to make a low-foaming solution	Cleaning of your evacuation equipment Only mix what you need so you eliminate waste and storage problems Cleaning of your eliminate	Box of 25 pouches Box of 60 pouches Box of 140 pouches	50036100— Unit Dose Pouch for Quarts Intro Pack— Box of 60 Pouches (Case Lot 24) 50036101— Unit Dose Pouch for Quarts Clinic Pack— Box of 140 Pouches (Case Lot 12) 50036102— Unit Dose Pouch for Gallons Clinic Pack— Box of 25 Pouches (Case Lot 24)

^{*} Each VK Solution Sold in Case Lots of Four (4) Gallons Only.